



## **TRUMAN TRAUMA X**

**Product Code:** TTR2000X

**Price:** Request a quote

**Size:** 79 x 45 x 28 cm

**Gross Weight:** 15kg

The TruMan Trauma X is ideal for practising a range of life threatening trauma procedures including; chest tube insertion, needle decompression, CPR, needle/surgical cricothyroidotomy and percutaneous tracheostomy. In addition to this, it also facilitates oral and nasal intubation training techniques.

### ***Ideal for Training in:***

- Chest tube insertion
- Needle decompression
- Cardio Pulmonary Resuscitation (CPR)
- Needle and surgical cricothyroidotomy
- Percutaneous tracheostomy
- Oral and nasal intubation
- Bag/mask ventilation
- Fibre optic intubation and video laryngoscopy
- Laryngeal mask airway and supraglottic device insertion
- Combi tube insertion
- Venus engorgement and tracheal deviation identification

### ***Features and Benefits:***

- Features the AirSim airway which provides realistic and durable training and comes with a 5 year warranty
- *Cricothyroidotomy and Tracheostomy*; A realistic larynx which features palpable cricoid landmarks, laryngeal cartilages and tracheal rings to provide positive user feedback and easy identification of the sternal notch
- *Chest tube insertion*: The manikin facilitates recognition of correct position, surgical incision, blunt dissection through chest wall, perforation of pleura, fluid simulation and finger sweep in the sixth and seventh intercostal space
- *Needle Decompression*: the user can create either right or left tension pneumothorax using a resuscitation bag; a hiss can be heard which the needle is successfully inserted into the second or fifth intercostal space

- *Cardio Pulmonary Resuscitation (CPR)*: easy identification of all anatomical landmarks, full head tilt with chin lift and jaw thrust to prepare the airway prior to resuscitation, lifelike recoil during compressions
- Cost effective replaceable consumables; realistic replaceable neck skin, replaceable larynx inserts which can facilitate 1 surgical cricothyroidotomy and 3 percutaneous tracheostomy procedures, needle decompression inserts which will facilitate up to 80+ needle incisions, chest drain inserts which have 3 intercostal spaces per insert
- A tongue bulb which can be inflated to create tongue edema providing a variance during training
- User feedback through visible chest rise
- Minimal maintenance and robust design which allows for efficient training
- A durable carry case for easy transportation and safe storage

***Package Contents:***

- 1 TruMan Trauma X (assembled and ready to use)
- 1 Black carrier case
- 1 USB pen drive user manual
- 1 Bottle of lubrication
- 2 Needle decompression inserts (as supplied on manikin)
- 2 Chest tube inserts ((as supplied on manikin)
- 3 Overlay neck skins (1 supplied on manikin and 2 spare)
- 5 larynx inserts (1 supplied on manikin and 4 spare)

