

US-10 US-10a US-10b



**NEW**

# Female Pelvic Ultrasound Phantom

**Practice both transvaginal and transabdominal ultrasound procedures in one**

Normal pregnancy unit is newly released with the pathological and ectopic pregnancy units!

**Product Supervision**

Charlotte Henningsen, MS, RT(R), RDMS, RVT, FSDMS, FAIUM  
 Chair & Professor - Sonography Department  
 Adventist University of Health Sciences



**FEATURES**

- 1 | Compatible with actual ultrasound probe
- 2 | Three types of interchangeable phantom inserts with different pathologies
- 3 | Realistic view of female external genital
- 4 | Excellent ultrasound image quality

**SKILLS**

- | Handling and manipulation of transvaginal and transabdominal transducers
- | Interpretation of sonogram
- | Visualization and localization of anatomies and pathologies

**ANATOMY** Three variations of training sets. Select the one which best meet your requirements.

US-10B: for Obstetrics	US-10A: Full set	US-10: for Gynecology
<p><b>Normal pregnancy unit (7 weeks)</b></p>	<p><b>Ectopic pregnancy unit</b></p>	<p><b>Pathological unit</b></p>

**DESCRIPTIONS**

**SPECIFICATIONS**

Torso size: W34 x D33 x H24cm  
 W13.5 x D13 x H9.5 inch

weight : 7.2kg/ 15.9 lbs

**SET INCLUDES**

US-10  
 Female Pelvic Ultrasound Phantom

- 1 | lower torso manikin
- 1 | ultrasound pathology unit
- 1 | ultrasound ectopic pregnancy unit
- 1 | storage case
- 1 | instruction manual

US-10A  
 Female Pelvic Ultrasound Phantom

- 1 | lower torso manikin
- 1 | ultrasound pathology unit
- 1 | ultrasound first trimester pregnancy unit
- 1 | ultrasound ectopic pregnancy unit
- 2 | storage case
- 1 | instruction manual

**MATERIALS**

Soft resin, Polyurethane elastomer  
 Latex free

US-10B  
 Obstetrics Pelvic Ultrasound Phantom

- 1 | lower torso manikin
- 1 | ultrasound first trimester pregnancy unit
- 1 | ultrasound ectopic pregnancy unit
- 1 | storage case
- 1 | instruction manual

**REPLACEMENT PARTS**

- 41909-030 | ultrasound pathology unit
- 41909-040 | ultrasound ectopic pregnancy unit
- 41909-050 | ultrasound first trimester pregnancy unit
- 41909-060 | skin for US-10/ 10A/ 10B

**RECOMMENDED DEVICES**

-  
 \* Ultrasound scanner is not included.