


Part No:

- 40201 Advanced
- 40200 Standard
- 40202 Basic

Breast Examination Trainer Advanced, Standard & Basic

**Learning Has Never
Been More Lifelike**



Pathologies can be placed in various predetermined location points & are easily changeable

Soft tissue looks & feels realistic

Hard torso available for bench top use

Dual purpose product:
bench top & hybrid use

Comfortable to wear for long periods of time during OSCEs

Easy for the user to fit & remove

Clavicular & axilla pads for lymph node placement

All components are latex free

Breast Examination Trainer - Advanced, Standard & Basic

This product provides a realistic platform for acquiring the delicate skills required for Clinical Breast Examination (CBE). Consisting of realistic soft tissue, with the versatility of interchangeable pathologies, it can be used either as a bench top teaching aid, or worn by a simulated patient during student assessment. An ideal product for any undergraduate programme running OSCEs, or healthcare professionals promoting best practice diagnostic skills to trainees & patients.

Designed in close collaboration with Mr Peter Donnelly of Torbay Hospital, UK. With the help & advice of staff at the Clinical Skills Resource Centre, University of Liverpool & the Clinical Practice Centre, St James's University Hospital, Leeds.



Advanced Breast Examination Trainer with Lymph Node Pads

Part No: 40201



Standard Breast Examination Trainer without Lymph Node Pads

Part No: 40200



Basic Breast Examination Trainer

Part No: 40202

Skills Gained:

- Clinical breast examination techniques (CBE & SBE)
- Professional-to-patient communication
- Identification of anatomical landmarks & lymph nodes (axillary, supra & infraclavicular)
- Location & diagnosis of pathologies



- Pathologies supplied: carcinomas: 2cm, 3cm, 5cm, cyst, fibrocystic disease, fibroadenoma
- Simulated carcinoma & fibroadenoma
- Benign cyst

+ Optional Extras:

- 40209 Breast Examination Carry Case

You'll Receive:

	Adv	Std	Basic
Wearable Examination Breasts	•	•	•
40204 Pathologies	•	•	•
40205 Breast Examination Inserts	•	•	•
40206 Breast Pathologies Supports	•	•	•
40207 Breast Back Plates	•	•	•
40208 Breast Examination Torso	•	•	-
40203 Lymph Node Pads + Lymph Nodes	•	-	-

Advanced

Packed weight (kg) & dimensions (mm)	5.5 kg	650 l	280 w	460 h
--------------------------------------	--------	-------	-------	-------

Standard

Packed weight (kg) & dimensions (mm)	5.3 kg	650 l	280 w	460 h
--------------------------------------	--------	-------	-------	-------

Basic

Packed weight (kg) & dimensions (mm)	1.5 kg	560 l	360 w	130 h
--------------------------------------	--------	-------	-------	-------

Patents granted

- US6336812 • EP099027 • DE69804793 • CA 2293585
- ZL98806331 • HK1023832 • 755575 Australia



Learn more at

limbsandthings.com

Limbs & Things Ltd. Sussex Street, St Philips, Bristol, BS2 0RA, UK

T: 00 44 (0)117 311 0500 F: 00 44 (0)117 311 0501

E: UK Customers - gbsales@limbsandthings.com

E: International Customers - sales@limbsandthings.com