



Clinical Male Pelvic Trainer (CMPT) Mk 2

Part No:

| 60951 Advanced

| 60950 Standard



Learn more at


limbsandthings.com

**Learning Has Never
Been More Lifelike**


Clinical Male Pelvic Trainer (CMPT) Mk 2
Advanced & Standard

The Clinical Male Pelvic Trainer Mk 2 presents key anatomical features, both externally and internally, for teaching ‘hands-on’ male pelvic examination and diagnosis. The uncluttered design of the model enables trainees to focus on essential anatomy and procedure. The Mk 2 model features new and improved pathologies which allow for a realistic and anatomically accurate training experience.


This product has been designed in close collaboration with the staff at the Clinical Skills Resource Centre, University of Liverpool; Imperial College, London; Guys Hospital, London; Southmead Hospital, Bristol; University of Southampton, UK.




Pelvis has realistic anatomy modelled from CT scan data




Skin surface is washable using soap and water




Soft tissue inserts are removable and replaceable




7 interchangeable modules available featuring pathologies of the genital area



Genitalia includes the penis (both circumcised and non-circumcised), Scrotum and testicles (includes vas deferens and epididymis)




Can be used in two positions (standing and supine)



Skills Gained


- Learning examination procedure
- Testicular examination
- Dry catheterisation (Foley catheter size 16 and retrograde catheters)
- Examination of normal anatomy
- Evaluate testicular abnormalities
- Evaluate abdominal / pelvic pain



CMPT Module 1:
Normal

Part No: **60952**


Packed weight (g) & dimensions (mm)				
720 g	600 l	210 w	160 h	



CMPT Module 2:
Varicocele

Part No: **60953**


Packed weight (g) & dimensions (mm)				
730 g	600 l	210 w	160 h	



CMPT Module 3:
Testicular Tumours and Penile Cancer

Part No: **60954**

Packed weight (g) & dimensions (mm)				
730 g	600 l	210 w	160 h	



CMPT Module 4:
Epididymal Cysts

Part No: **60955**


Packed weight (g) & dimensions (mm)				
740 g	600 l	210 w	160 h	



CMPT Module 5:
Hydrocele

Part No: **60956**


Packed weight (g) & dimensions (mm)				
835 g	600 l	210 w	160 h	



CMPT Module 6:
Orchitis/ Epididymo-orchitis

Part No: **60957**

Packed weight (g) & dimensions (mm)				
810 g	600 l	210 w	160 h	



CMPT Module 7:
Indirect Inguinal Hernia

Part No: **60958**

Packed weight (g) & dimensions (mm)				
775 g	600 l	210 w	160 h	

Clinical Male Pelvic Trainer (CMPT) Mk 2

Part No:

60951 Advanced
60950 Standard



The genitalia includes the penis (both circumcised and non-circumcised), scrotum, testicles (includes vas deferens and epididymis)

Anatomy of penis:

- glans and shaft normal

Anatomy of groin with clear anatomical landmarks:

- ASIS (anterior superior iliac spine)
- pubic symphysis and tubercles

Anatomy of scrotum area:

- scrotum
- testicles with epididymis and vas deferens

Pathologies of genital area (Advanced model only):

- varicocele (Module 2)
- testicular tumour & penile cancer (Module 3)
- epididymal cysts (Module 4)
- hydrocele (Module 5)
- orchitis/epididymo-orchitis (Module 6)
- indirect Inguinal hernia (Module 7)



Skills Gained:

- Learning examination procedure
- Testicular examination
- Dry catheterisation (Foley catheter size 16 and retrograde catheters)
- Examination of normal anatomy
- Evaluate testicular abnormalities
- Evaluate abdominal / pelvic pain

Product Qualities



- Soft tissue inserts are removable and replaceable:
 - genitals (penis, perineum, & testicles)
 - abdominal wall



- Can be used in two positions (standing and supine)



- Skin surface is washable using soap and water

Advanced - 60951

You'll Receive:

- **60952** 1 CMPT Module 1: Normal
- **60953** 1 CMPT Module 2: Varicocele
- **60954** 1 CMPT Module 3: Tumours
- **60955** 1 CMPT Module 4: Epididymal Cysts
- **60956** 1 CMPT Module 5: Hydrocele
- **60957** 1 CMPT Module 6: Orchitis/Epididymo-orchitis
- **60958** 1 CMPT Module 7: Indirect Inguinal Hernia
- **60959** 1 CMPT Abdominal Insert
- **60960** 2 CMPT Pubic Bone
- **60961** 2 CMPT Foreskins Pack (each pack contains 2 foreskins and 3 lubricant sachets)
- 1 Lower Torso
- 1 Support Base
- 1 Carry Case
- 1 LED Pen Torch

Trainer

Packed weight (kg) & dimensions (mm) 13.2 kg 650 l 290 w 480 h

Penis

Packed weight (kg) & dimensions (mm) 6.9 kg 550 l 290 w 480 h

Standard - 60950

You'll Receive:

- **60952** 1 CMPT Module 1 Normal
- **60959** 1 CMPT Abdominal Insert
- **60960** 1 CMPT Pubic Bone
- **60961** 1 CMPT Foreskins Pack (contains 2 foreskins and 3 lubricant sachets)
- 1 Lower Torso
- 1 Support Base
- 1 Carry Case

Packed weight (kg) & dimensions (mm) 13.2 kg 650 l 290 w 480 h

Patents granted

- US6336812 • EP099027 • DE69804793 • CA 2293585
- ZL98806331 • HK1023832 • 755575 Australia



Learn more at:

limbsandthings.com

Limbs & Things Ltd. Sussex Street, St Philips, Bristol, BS2 0RA, UK

T: 00 44 (0)117 311 0500 F: 00 44 (0)117 311 0501

E: UK Customers - gbsales@limbsandthings.com

E: International Customers - sales@limbsandthings.com