




Fundamentals of Laparoscopic Surgery™ Trainer System & Accessories



*Validated
Laparoscopic Skills Training*

Fundamentals of Laparoscopic Surgery™ Trainer System

Working closely with the SAGES-ACS FLS committee, we have revolutionised the FLS trainer system. Introducing an updated series of laparoscopic training products exclusively designed and endorsed for the FLS Program.


The systems come complete with a training box and TV camera alongside a range of tasks designed to develop the skills necessary to perform the fundamentals of laparoscopic surgery.

FLS Trainer Box with TV Camera

This product is economical and portable, making it ideal for training centres.

Note that users must provide their own TV screen.

Product No: **50302**



FLS Trainer System with TV Camera

The FLS Trainer System incorporates all of the features of the FLS Trainer Box but also includes a 19" Television Screen and a Mobile Sim Cart. This product offers greater flexibility and is extremely cost-effective option for training centres.



+



+



Product No: 50304

FLS All-In-One Trainer System with TV Camera

The FLS All-In-One Trainer System incorporates all of the features of the FLS Trainer System and also includes ligating loops, sutures and an FLS Approved Instrument Kit.

This product provides a training centre with everything required to run an FLS Training program.



+



+



Product No: 50306

Check our website for more information about Optional Consumables & Instruments



For further information
on the FLS training system visit

fls-products.com

*Learning
has never
been more
lifelike*

FLS Trainer System

	FLS Trainer Box with TV Camera	FLS Trainer System with TV Camera	FLS All-In-One Trainer System
	<i>Product No:</i> 50302	<i>Product No:</i> 50304	<i>Product No:</i> 50306
Trainer Box	✓	✓	✓
Camera - TV (with lead)	✓	✓	✓
Light Strip	✓	✓	✓
Power Cable	✓	✓	✓
Simulated Skin	✓	✓	✓
2 x String & Alligator Clips	✓	✓	✓
2 x Trocar	✓	✓	✓
1 x Peg Board with 6 Triangles	✓	✓	✓
1 x Jumbo Clip Retainer	✓	✓	✓
25 x Single Circle Gauze	✓	✓	✓
25 x Double Circle Gauze	✓	✓	✓
3 x Foam Organs	✓	✓	✓
50 x Penrose Drains	✓	✓	✓
3 x Single Use Ligation Loops	•	•	✓
5 x Multi-Use Ligation Loops	•	•	✓
50 x Sutures 2/0 Silk 90cm, 26mm Needle	•	•	✓
19" TV	•	✓	✓
Mobile SimCart	•	✓	✓
1 x FLS Approved Instrument Kit	•	•	✓



limbsandthings.com

Limbs & Things Ltd. Sussex Street, St Philips, Bristol, BS2 0RA, UK

T: 00 44 (0)117 311 0500 F: 00 44 (0)117 311 0501 E: gbsales@limbsandthings.com

UK September 2016