



Uterine Involution Assessment Simulator

Practice highly realistic uterine involution assessments via a soft, durable and elastic postpartum torso





Palpation of uterine fundus




Measurement of fundal height

FEATURES

- 1 | Soft and elastic abdominal wall that permits highly realistic palpation to assess the degree of uterine involution and to measure fundal height.
- 2 | Four interchangeable uterine inserts

SKILLS

- Assessment of uterus in early puerperium
- Fundal height measurement
- Palpation of uterine fundus
- Perineal cleansing and assessment
- Postpartum massage
- Assessment of anal prolapse



SPECIFICATIONS

Size: W28 x D42 x H26 cm W11 x D16.5 x 10.2 cm Weight: 7.4kg/ 16.3 lbs

SET INCLUDES

- 1 puerperal lower torso
- set of uterine inserts (4 variations)
- 1 talcum powder
- instruction manual

DESCRIPTIONS

MATERIALS

Soft resin Latex free

RECOMMENDED DEVICES

IREPLACEMENT PARTS

11418-010 1 skin cover (for uterine involution assessment)

11418-020 1 set of uterine inserts (4 variations)

