



**Comprehensive training set with the platform of the maternal body for skills training in pelvic vaginal examination during labor and delivery assistance.**

### Features

#### For all training skills

- Genital area made of innovative material that excels in elasticity, durability, and resilience, offers repetitive training.



#### Delivery assistance

- Anatomically correct vaginal canal, ischial spines, pubic symphysis and coccyx that facilitate training with anatomical understanding.
- Elastic external genitalia, realistic shape of vaginal canal and the holding handle for the maternal body to facilitate simulating a smooth delivery.



#### Cervical examination

- Five interchangeable stages of cervical dilation module can be installed with a single action. Descent of the fetal head can also be set.



- Anatomically correct landmarks including the ischial spine and the pubic symphysis
- Anterior fontanelle can be touched to determine the fetal position.

#### Perineorrhaphy (perineal repair)

- Repetitive suture training for the first-degree perineal tears.



\* Suturing instruments are not included.

## Cervical examination

### Training Skills

- Cervical examination
  - ▶ Determination of station using landmarks of the ischial spine and the pubic joint to assess the descent of the body
  - ▶ Assessment of cervical dilation
- Speculum examination
- Procedures of observation and cervical cytology
- Procedures of urinary catheterization



Pelvic vaginal examination



Urinary catheterization

\* For replacement parts, refer to page 8.

## Delivery assistance

### Training Skills

- Delivery assistance for different delivery positions (dorsal, all-fours)
- Perineal protection
- Delivery of fetus (forceps delivery, vacuum extraction)
- Clamping, tying and cutting of the umbilical cord
- Delivery of the placenta
- Insertion of a urinary catheter



Perineal protection



Delivery assistance



Cutting the umbilical cord



Delivery of the placenta

\* For replacement parts, refer to page 9.

## Perineorrhaphy (perineal repair)

### Training Skills

- Perineorrhaphy (perineal repair)
  - ▶ Single interrupted suture
  - ▶ Vertical mattress suture



### Specifications :

Common Components:	Parts for pelvic examination	Parts for Delivery Assistance	Parts for Perineorrhaphy
1 maternal body	1 vagina unit	1 fetus	1 genitalia unit
1 genitalia unit (for pelvic examination and delivery assistance)	1 set of cervical dilation inserts (5 variations)	1 placenta	(for perineorrhaphy)
1 talcum powder	1 holder for the cervical dilation inserts	3 umbilical cords for omphalotomy	
3 lubricant (150ml)		1 tray	
1 instruction manual		5 velamen sheets	